

For advanced video laryngoscopy, the choice has never been more clear.

The GlideScope AVL Reusable features a digital color monitor and digital camera, providing a consistently clear, DVD-quality view, enabling swift intubation.

The AVL system gives you a choice of single-use and reusable options designed for a wide range of patient weights and sizes—ideal for demanding OR, ED, and ICU settings.

Add real-time recording and an onboard video tutorial, and these versatile products are a clear choice in advanced video laryngoscopy.

#### AVL=ADVANCED AIRWAY VIEWS


- ▶ Digital color monitor with DVD-quality resolution
- ▶ Reusable video laryngoscope with digital-resolution camera
- ▶ Advanced resolution output for OR theater viewing
- ▶ Anti-fogging mechanism resists lens clouding and secretions

#### AVL=ADVANCED USER FEATURES

- ▶ Integrated, real-time recording helps confirm tube placement and facilitate teaching
- ▶ Onboard *4-Step Technique* tutorial makes AVL easy to use, learn, and teach
- ▶ Ergonomic design and rugged construction for multiple clinical applications
- ▶ Reusable blades in 4 sizes: GVL® 2, 3, 4, 5



## Advanced Video Laryngoscopy: Real-time recording. DVD clarity.



Illustrated airway view.

Applications

- First-use intubations, replacing direct laryngoscopy
- Normal or restricted oropharyngeal views
- Cormack-Lehane Grades I–IV laryngeal views
- Trauma airways—excellent when dealing with blood and secretions in the airway
- Nasal tracheal intubation
- Video-guided tube exchange in the ICU

System Components

- Digital, DVD-quality color monitor
- Reusable video laryngoscope with digital camera
- Reusable AVL blades in a choice of four (GVL®) sizes
- GlideRite® Rigid Stylet
- Choice of mobile stand or IV pole mount
- In-Service CD including User’s Manual and Quick Reference Cards

GlideScope Products and Systems: Designed for 1<sup>st</sup> Pass Success

Reusable Size	GVL 2	GVL 3	GVL 4	GVL 5
Patient Weight	4 - 20 kg	10 kg - Adult	40 kg - Morbidly Obese	40 kg - Morbidly Obese
GlideRite Rigid Stylet	For use with endotracheal tubes 6.0 mm and larger.			

GVL 2	GVL 3	GVL 4	GVL 5*
			
Length (blade tip to handle): 51 mm Height at camera: 11 mm Width at camera: 14 mm	Length (blade tip to handle): 82 mm Height at camera: 16 mm Width at camera: 20 mm	Length (blade tip to handle): 103 mm Height at camera: 14 mm Width at camera: 27 mm	Length (blade tip to handle): 103 mm Height at camera: 14 mm Width at camera: 27 mm

Note: System configurations vary.

\*The GVL 5 is designed to accommodate anatomical anomalies sometimes associated with bariatric patients.

AVL PRODUCT CONFIGURATION	PRODUCT OPTIONS
 <p><b>On Mobile Stand</b> A stand-alone configuration for use in a variety of hospital settings</p> <p>Power source: battery and/or AC power</p> <p><b>Also available in IV Pole Mount</b></p>	 <p><b>GlideScope Direct</b> The GlideScope Direct intubation trainer is designed to facilitate instruction of classic direct laryngoscopy (DL).</p> <p><b>Available in Single-use and Reusable</b></p>



GlideScope, the GlideScope symbol, GVL, GlideRite, Verathon, and the Verathon Torch symbol are trademarks of Verathon Inc. © 2011 Verathon Inc.

GlideScope® video laryngoscope systems are CE marked in accordance with the Medical Device Directive, and the Verathon Inc. quality system is Quality System Certified to ISO 13485:2003 standards.

EUSM-0367-00-86

**Verathon Medical Europe B.V.**  
Linnaeusweg 11 / 3401 MS IJsselstein / The Netherlands  
Tel. +31 (0)30 68 70 570 / Fax +31 (0)30 68 70 512  
verathon.eu